

First Stage: Natrelle® 133S Smooth Tissue Expanders

Select an expander below based on the width measurement and find the complementary range of *Natrelle* INSPIRA® Smooth Round Gel Breast Implants to the right in the same color.*

Style FV-T
Full Height, Variable Projection



133S
(Smooth)

WIDTH (CM)	HEIGHT (CM)	PROJECTION (CM)	FILL VOLUME (CC)	CATALOG NUMBER
11.0	11.5	5.0	300	133S-FV-11-T
12.0	12.5	5.3	400	133S-FV-12-T
13.0	13.5	5.7	500	133S-FV-13-T
14.0	14.5	6.2	600	133S-FV-14-T
15.0	15.5	6.7	750	133S-FV-15-T
16.0	16.5	6.8	850	133S-FV-16-T

Style MV-T
Moderate Height, Variable Projection



133S
(Smooth)

WIDTH (CM)	HEIGHT (CM)	PROJECTION (CM)	FILL VOLUME (CC)	CATALOG NUMBER
11.0	10.0	4.9	250	133S-MV-11-T
12.0	11.0	5.2	300	133S-MV-12-T
13.0	12.0	5.6	400	133S-MV-13-T
14.0	13.0	6.0	500	133S-MV-14-T
15.0	14.0	6.3	600	133S-MV-15-T
16.0	15.0	6.6	700	133S-MV-16-T

Style SV-T
Short Height, Variable Projection



133S
(Smooth)

WIDTH (CM)	HEIGHT (CM)	PROJECTION (CM)	FILL VOLUME (CC)	CATALOG NUMBER
11.0	7.6	4.7	200	133S-SV-11-T
12.0	8.4	5.3	250	133S-SV-12-T
13.0	9.1	5.7	300	133S-SV-13-T
14.0	9.7	6.0	375	133S-SV-14-T
15.0	10.5	6.3	450	133S-SV-15-T
16.0	11.3	6.3	550	133S-SV-16-T

Note: *Natrelle*® 133S Smooth Tissue Expanders have 6 suture tabs (except for Style SV, which has 4).
*For implant widths from 11.0 cm to 16.5 cm.

Style FX-T
Full Height, Extra Projection



133S
(Smooth)

WIDTH (CM)	HEIGHT (CM)	PROJECTION (CM)	FILL VOLUME (CC)	CATALOG NUMBER
11.0	11.5	5.9	350	133S-FX-11-T
12.0	12.5	6.3	450	133S-FX-12-T
13.0	13.5	6.7	550	133S-FX-13-T
14.0	14.5	7.1	650	133S-FX-14-T
15.0	15.5	7.6	800	133S-FX-15-T
15.5	16.0	8.1	950	133S-FX-16-T

Style MX-T
Moderate Height, Extra Projection



133S
(Smooth)

WIDTH (CM)	HEIGHT (CM)	PROJECTION (CM)	FILL VOLUME (CC)	CATALOG NUMBER
11.0	10.0	5.9	300	133S-MX-11-T
12.0	11.0	6.3	400	133S-MX-12-T
13.0	12.0	6.7	500	133S-MX-13-T
14.0	13.0	7.1	600	133S-MX-14-T
15.0	14.0	7.6	700	133S-MX-15-T
15.5	14.5	8.1	850	133S-MX-16-T

Style SX-T
Short Height, Extra Projection



133S
(Smooth)

WIDTH (CM)	HEIGHT (CM)	PROJECTION (CM)	FILL VOLUME (CC)	CATALOG NUMBER
11.0	9.0	5.9	250	133S-SX-11-T
12.0	10.0	6.3	350	133S-SX-12-T
13.0	11.0	6.7	400	133S-SX-13-T
14.0	12.0	7.1	500	133S-SX-14-T
15.0	13.0	7.6	650	133S-SX-15-T
15.5	13.5	8.1	800	133S-SX-16-T



Second Stage: Natrelle INSPIRA® SoftTouch Smooth Round Gel Breast Implants



These color-coded charts are intended as guides to assist you in selecting complementary expanders and implants.

Natrelle INSPIRA® Style SSL
Smooth, Low Profile



WIDTH (CM)	PROJECTION (CM)	VOLUME (CC)	Smooth SoftTouch (SSL)
10.00	2.0	110	SSL-110
10.25	2.0	125	SSL-125
10.50	2.1	140	SSL-140
11.00	2.2	170	SSL-170
11.50	2.3	200	SSL-200
12.00	2.4	230	SSL-230
12.50	2.5	260	SSL-260
13.00	2.6	290	SSL-290
13.50	2.7	320	SSL-320
14.00	2.8	350	SSL-350
14.50	2.9	380	SSL-380
15.00	3.0	410	SSL-410
15.50	3.1	460	SSL-460
16.00	3.2	510	SSL-510
16.50	3.3	560	SSL-560
17.00	3.4	610	SSL-610

Natrelle INSPIRA® Style SSLP
Smooth, Low-Plus Profile



WIDTH (CM)	PROJECTION (CM)	VOLUME (CC)	Smooth SoftTouch (SSLP)
9.50	2.5	125	SSLP-125
10.00	2.6	145	SSLP-145
10.25	2.7	165	SSLP-165
11.00	2.8	190	SSLP-190
11.25	2.9	205	SSLP-205
11.50	3.0	220	SSLP-220
11.75	3.1	235	SSLP-235
12.00	3.1	250	SSLP-250
12.25	3.2	265	SSLP-265
12.50	3.3	280	SSLP-280
12.75	3.3	300	SSLP-300
13.00	3.4	320	SSLP-320
13.25	3.5	340	SSLP-340
13.50	3.5	360	SSLP-360
14.00	3.6	400	SSLP-400
14.50	3.8	440	SSLP-440
15.00	3.9	490	SSLP-490
15.50	4.0	540	SSLP-540
16.00	4.2	590	SSLP-590
16.50	4.3	640	SSLP-640

Natrelle INSPIRA® Style SSM
Smooth, Moderate Profile



WIDTH (CM)	PROJECTION (CM)	VOLUME (CC)	Smooth SoftTouch (SSM)
9.50	3.0	140	SSM-140
10.00	3.3	175	SSM-175
10.25	3.4	195	SSM-195
10.50	3.5	210	SSM-210
11.00	3.7	240	SSM-240
11.25	3.8	255	SSM-255
11.50	3.9	275	SSM-275
11.75	3.9	295	SSM-295
12.00	4.0	310	SSM-310
12.25	4.1	330	SSM-330
12.50	4.2	345	SSM-345
12.75	4.2	360	SSM-360
13.00	4.3	375	SSM-375
13.25	4.4	405	SSM-405
13.50	4.5	445	SSM-445
14.00	4.6	485	SSM-485
14.50	4.8	520	SSM-520
14.75	4.9	560	SSM-560
15.00	4.9	600	SSM-600
15.50	5.0	640	SSM-640
16.00	5.2	685	SSM-685
16.50	5.0	755	SSM-755

Natrelle INSPIRA® Style SSF
Smooth, Full Profile



WIDTH (CM)	PROJECTION (CM)	VOLUME (CC)	Smooth SoftTouch (SSF)
9.50	4.0	180	SSF-180
9.75	4.1	200	SSF-200
10.00	4.2	220	SSF-220
10.50	4.3	240	SSF-240
11.00	4.4	265	SSF-265
11.25	4.5	295	SSF-295
11.40	4.6	325	SSF-325
11.50	4.7	335	SSF-335
11.75	4.8	345	SSF-345
12.00	4.9	365	SSF-365
12.25	5.0	385	SSF-385
12.50	5.1	415	SSF-415
12.75	5.3	450	SSF-450
13.00	5.4	485	SSF-485
13.25	5.5	520	SSF-520
13.50	5.7	560	SSF-560
14.00	5.8	605	SSF-605
14.50	5.9	650	SSF-650
14.75	6.0	695	SSF-695
15.00	6.1	745	SSF-745
15.50	6.0	770	SSF-770

Natrelle INSPIRA® Style SSX
Smooth, Extra-Full Profile



WIDTH (CM)	PROJECTION (CM)	VOLUME (CC)	Smooth SoftTouch (SSX)
9.50	4.2	200	SSX-200
10.00	4.5	225	SSX-225
10.25	4.8	255	SSX-255
10.50	5.0	285	SSX-285
10.75	5.2	310	SSX-310
11.00	5.4	340	SSX-340
11.25	5.6	375	SSX-375
11.50	5.8	400	SSX-400
11.75	5.9	420	SSX-420
12.00	6.0	445	SSX-445
12.25	6.1	470	SSX-470
12.50	6.1	495	SSX-495
12.75	6.2	525	SSX-525
13.00	6.2	545	SSX-545
13.25	6.3	560	SSX-560
13.50	6.4	580	SSX-580
13.75	6.4	615	SSX-615
14.00	6.5	650	SSX-650
14.25	6.6	700	SSX-700
14.50	6.7	750	SSX-750
15.00	6.7	800	SSX-800

Natrelle® Breast Implants

IMPORTANT SAFETY INFORMATION (continued)

INDICATIONS

Natrelle® Breast Implants are indicated for women for the following:

- **Breast augmentation for women at least 22 years old for silicone-filled implants and breast augmentation for women at least 18 years old for saline-filled implants.** This includes primary breast augmentation to increase the breast size, as well as revision surgery to correct or improve the result of a primary breast augmentation surgery
- **Breast reconstruction.** This includes primary reconstruction to replace breast tissue that has been removed due to cancer or trauma or that has failed to develop properly due to a severe breast abnormality. Breast reconstruction also includes revision surgery to correct or improve the result of a primary breast reconstruction surgery

CONTRAINDICATIONS

Breast implant surgery should not be performed in:

- Women with active infection anywhere in their body
- Women with existing cancer or precancer of their breast who have not received adequate treatment for those conditions
- Women who are currently pregnant or nursing

ADDITIONAL WARNINGS

- See Boxed Warning
- **Avoid damage during surgery:** Care should be taken to avoid the use of excessive force and to minimize handling of the implant. Forcing of implants through too small an opening or applying concentrated localized pressure on the implants may result in localized weakening of the breast implant shell, potentially leading to shell damage and possible implant rupture. An incision should be of appropriate length to accommodate the style, size, and profile of the implants. Use care when using surgical instruments in proximity with the breast implant
- Follow recommended fill volumes for saline implants to decrease possibility of shell wrinkling and crease-fold failure

PRECAUTIONS

Safety and effectiveness have not been established in patients with the following:

- Autoimmune diseases (eg, lupus and scleroderma)
- A compromised immune system (eg, currently receiving immunosuppressive therapy)
- Planned chemotherapy or radiation following breast implant placement
- Conditions or medications that interfere with wound healing and blood clotting
- Reduced blood supply to breast tissue
- Clinical diagnosis of depression or other mental health disorders, including body dysmorphic disorder and eating disorders. Please discuss any history of mental health disorders prior to surgery. Patients with a diagnosis of depression, or other mental health disorders, should wait until resolution or stabilization of these conditions prior to undergoing breast implantation surgery

ADVERSE EVENTS

Possible adverse events with breast implant surgery include implant rupture with silicone implants, implant deflation with saline-filled implants, capsular contracture, reoperation, implant removal, pain, changes in nipple and breast sensation, infection, scarring, asymmetry, wrinkling, implant displacement/migration, implant palpability/visibility, breastfeeding complications, hematoma/seroma, implant extrusion, necrosis, delayed wound healing, infection, breast tissue atrophy/chest wall deformity, calcium deposits, and lymphadenopathy. Other systemic conditions have been reported with breast implants.

For more information, please see the full Directions for Use at www.allergan.com/products.

To report a problem with *Natrelle*® Breast Implants, please call Allergan® at 1-800-624-4261.

The sale and distribution of this device is restricted to users and/or user facilities that provide information to patients about the risks and benefits of this device in the form and manner specified in the approved labeling provided by Allergan®.

Matchup Guide

Match *Natrelle*® 133S Smooth Tissue Expanders with *Natrelle* INSPIRA® Cohesive and SoftTouch Round Gel Breast Implants.

Natrelle INSPIRA® Smooth Round Gel Breast Implants



Natrelle® 133S Smooth Tissue Expander



FOURTE™ Expander Fill System
Included with all *Natrelle*® 133S Smooth Tissue Expanders



Please see Important Safety Information inside, including Boxed Warning.

First Stage: Natrelle® 133S Smooth Tissue Expanders

Select an expander below based on the width measurement and find the complementary range of Natrelle INSPIRA® Smooth Round Gel Breast Implants to the right in the same color.*

Style FV-T Full Height, Variable Projection					133S (Smooth)
WIDTH (CM)	HEIGHT (CM)	PROJECTION (CM)	FILL VOLUME (CC)	CATALOG NUMBER	
11.0	11.5	5.0	300	133S-FV-11-T	
12.0	12.5	5.3	400	133S-FV-12-T	
13.0	13.5	5.7	500	133S-FV-13-T	
14.0	14.5	6.2	600	133S-FV-14-T	
15.0	15.5	6.7	750	133S-FV-15-T	
16.0	16.5	6.8	850	133S-FV-16-T	

Style MV-T Moderate Height, Variable Projection					133S (Smooth)
WIDTH (CM)	HEIGHT (CM)	PROJECTION (CM)	FILL VOLUME (CC)	CATALOG NUMBER	
11.0	10.0	4.9	250	133S-MV-11-T	
12.0	11.0	5.2	300	133S-MV-12-T	
13.0	12.0	5.6	400	133S-MV-13-T	
14.0	13.0	6.0	500	133S-MV-14-T	
15.0	14.0	6.3	600	133S-MV-15-T	
16.0	15.0	6.6	700	133S-MV-16-T	

Style SV-T Short Height, Variable Projection					133S (Smooth)
WIDTH (CM)	HEIGHT (CM)	PROJECTION (CM)	FILL VOLUME (CC)	CATALOG NUMBER	
11.0	7.6	4.7	200	133S-SV-11-T	
12.0	8.4	5.3	250	133S-SV-12-T	
13.0	9.1	5.7	300	133S-SV-13-T	
14.0	9.7	6.0	375	133S-SV-14-T	
15.0	10.5	6.3	450	133S-SV-15-T	
16.0	11.3	6.3	550	133S-SV-16-T	

Note: Natrelle® 133S Smooth Tissue Expanders have 6 suture tabs (except for Style SV, which has 4).

*For implant widths from 11.0 cm to 16.5 cm.

Style FX-T Full Height, Extra Projection					133S (Smooth)
WIDTH (CM)	HEIGHT (CM)	PROJECTION (CM)	FILL VOLUME (CC)	CATALOG NUMBER	
11.0	11.5	5.9	350	133S-FX-11-T	
12.0	12.5	6.3	450	133S-FX-12-T	
13.0	13.5	6.7	550	133S-FX-13-T	
14.0	14.5	7.1	650	133S-FX-14-T	
15.0	15.5	7.6	800	133S-FX-15-T	
15.5	16.0	8.1	950	133S-FX-16-T	

Style MX-T Moderate Height, Extra Projection					133S (Smooth)
WIDTH (CM)	HEIGHT (CM)	PROJECTION (CM)	FILL VOLUME (CC)	CATALOG NUMBER	
11.0	10.0	5.9	300	133S-MX-11-T	
12.0	11.0	6.3	400	133S-MX-12-T	
13.0	12.0	6.7	500	133S-MX-13-T	
14.0	13.0	7.1	600	133S-MX-14-T	
15.0	14.0	7.6	700	133S-MX-15-T	
15.5	14.5	8.1	850	133S-MX-16-T	

Style SX-T Short Height, Extra Projection					133S (Smooth)
WIDTH (CM)	HEIGHT (CM)	PROJECTION (CM)	FILL VOLUME (CC)	CATALOG NUMBER	
11.0	9.0	5.9	250	133S-SX-11-T	
12.0	10.0	6.3	350	133S-SX-12-T	
13.0	11.0	6.7	400	133S-SX-13-T	
14.0	12.0	7.1	500	133S-SX-14-T	
15.0	13.0	7.6	650	133S-SX-15-T	
15.5	13.5	8.1	800	133S-SX-16-T	

Natrelle®
133S smooth tissue expander

Second Stage: Natrelle INSPIRA® Cohesive Smooth Round Gel Breast Implants

These color-coded charts are intended as guides to assist you in selecting complementary expanders and implants.

Natrelle INSPIRA® Style SCL Smooth, Low Profile			
WIDTH (CM)	PROJECTION (CM)	VOLUME (CC)	Smooth Cohesive (SCL)
---	---	---	---
10.00	2.0	110	SCL-110
10.25	2.0	125	SCL-125
10.50	2.1	140	SCL-140
---	---	---	---
11.00	2.2	170	SCL-170
---	---	---	---
11.50	2.3	200	SCL-200
---	---	---	---
12.00	2.4	230	SCL-230
---	---	---	---
12.50	2.5	260	SCL-260
---	---	---	---
13.00	2.6	290	SCL-290
---	---	---	---
13.50	2.7	320	SCL-320
---	---	---	---
14.00	2.8	350	SCL-350
---	---	---	---
14.50	2.9	380	SCL-380
---	---	---	---
15.00	3.0	410	SCL-410
15.50	3.1	460	SCL-460
16.00	3.2	510	SCL-510
16.50	3.3	560	SCL-560
17.00	3.4	610	SCL-610

Natrelle® 133S Smooth Tissue Expanders With MAGNA-SITE® Injection Sites **IMPORTANT SAFETY INFORMATION**

INDICATIONS

- Natrelle® 133S Smooth Tissue Expanders are indicated for:
- Breast reconstruction following mastectomy
 - Treatment of underdeveloped breasts
 - Treatment of soft tissue deformities

IMPORTANT SAFETY INFORMATION CONTRAINDICATIONS

Natrelle® 133S Smooth Tissue Expanders should not be used in patients:

- Who already have implanted devices that would be affected by a magnetic field (eg, pacemakers, drug infusion devices, and artificial sensing devices)
- Whose tissue at the expansion site is determined to be unsuitable
- Who have an active infection or a residual gross tumor at the expansion site
- Who are undergoing adjuvant radiation therapy

Natrelle INSPIRA® Style SCLP Smooth, Low-Plus Profile			
WIDTH (CM)	PROJECTION (CM)	VOLUME (CC)	Smooth Cohesive (SCLP)
9.50	2.5	125	SCLP-125
10.00	2.6	145	SCLP-145
---	---	---	---
10.50	2.7	165	SCLP-165
---	---	---	---
11.00	2.8	190	SCLP-190
11.25	2.9	205	SCLP-205
---	---	---	---
11.50	3.0	220	SCLP-220
11.75	3.1	235	SCLP-235
12.00	3.1	250	SCLP-250
12.25	3.2	265	SCLP-265
12.50	3.3	280	SCLP-280
12.75	3.3	300	SCLP-300
13.00	3.4	320	SCLP-320
13.25	3.5	340	SCLP-340
13.50	3.5	360	SCLP-360
---	---	---	---
14.00	3.6	400	SCLP-400
---	---	---	---
14.50	3.8	440	SCLP-440
---	---	---	---
15.00	3.9	490	SCLP-490
15.50	4.0	540	SCLP-540
16.00	4.2	590	SCLP-590
16.50	4.3	640	SCLP-640
---	---	---	---

- Whose physiological condition (eg, sensitive over- or underlying anatomy, obesity, smoking, diabetes, autoimmune disease, hypertension, chronic lung or severe cardiovascular disease, or osteogenesis imperfecta) or use of certain drugs (including those that interfere with blood clotting or affect tissue viability) poses an unduly high risk of surgical and/or postoperative complications
- Who are psychologically unsuitable

WARNINGS

- Do not use Natrelle® 133S Smooth Tissue Expanders in patients who already have implanted devices that would be affected by a magnetic field (see *Contraindications*) because the MAGNA-SITE® integrated injection site contains a strong rare-earth, permanent magnet. Diagnostic testing with Magnetic Resonance Imaging (MRI) is contraindicated in patients with Natrelle® 133S Smooth Tissue Expanders in place
- Do not alter the tissue expander or use adulterated fill. Fill only with sterile saline for injection as described in INSTRUCTIONS FOR USE. Do not expose to contaminants
- Do not expand if the pressure will compromise wound healing or vasculature of overlying tissue, or beyond patient

Natrelle INSPIRA® Style SCM Smooth, Moderate Profile			
WIDTH (CM)	PROJECTION (CM)	VOLUME (CC)	Smooth Cohesive (SCM)
9.50	3.0	140	SCM-140
10.00	3.3	175	SCM-175
10.25	3.4	195	SCM-195
10.50	3.5	210	SCM-210
---	---	---	---
11.00	3.7	240	SCM-240
11.25	3.8	255	SCM-255
---	---	---	---
11.50	3.9	275	SCM-275
11.75	3.9	295	SCM-295
12.00	4.0	310	SCM-310
12.25	4.1	330	SCM-330
12.50	4.2	345	SCM-345
12.75	4.2	360	SCM-360
13.00	4.3	375	SCM-375
13.25	4.4	405	SCM-405
13.50	4.5	445	SCM-445
---	---	---	---
14.00	4.6	485	SCM-485
---	---	---	---
14.50	4.8	520	SCM-520
14.75	4.9	560	SCM-560
15.00	4.9	600	SCM-600
15.50	5.0	640	SCM-640
16.00	5.2	685	SCM-685
16.50	5.0	755	SCM-755
---	---	---	---

- tissue tolerance. Stop filling immediately if tissue damage, wound dehiscence, abnormal skin pallor, erythema, edema, pain, or tenderness are observed
- Do not reuse explanted products
- Active infection anywhere may increase risk of periprosthetic infection. Do not expose the tissue expander or injection needles to contaminants. Postoperative infections should be treated aggressively
- Adverse reactions may require premature explantation
- When using suturing tabs be careful to avoid piercing the shell. Use a new one if damage occurs
- Natrelle® 133S Smooth Tissue Expanders are temporary devices and are not to be used for permanent implantation or beyond 6 months. Tissue expansion in breast reconstruction typically requires 4 to 6 months

Natrelle INSPIRA® Style SCF Smooth, Full Profile			
WIDTH (CM)	PROJECTION (CM)	VOLUME (CC)	Smooth Cohesive (SCF)
9.50	4.0	180	SCF-180
9.75	4.1	200	SCF-200
10.00	4.2	220	SCF-220
10.50	4.3	240	SCF-240
---	---	---	---
11.00	4.4	265	SCF-265
11.25	4.5	295	SCF-295
11.40	4.6	325	SCF-325
11.50	4.7	335	SCF-335
11.75	4.8	345	SCF-345
12.00	4.9	365	SCF-365
12.25	5.0	385	SCF-385
12.50	5.1	415	SCF-415
12.75	5.3	450	SCF-450
13.00	5.4	485	SCF-485
13.25	5.5	520	SCF-520
13.50	5.7	560	SCF-560
---	---	---	---
14.00	5.8	605	SCF-605
---	---	---	---
14.50	5.9	650	SCF-650
14.75	6.0	695	SCF-695
15.00	6.1	745	SCF-745
15.50	6.0	770	SCF-770
---	---	---	---
---	---	---	---
---	---	---	---
---	---	---	---

PRECAUTIONS

Active infections may need to be treated and resolved before surgery. Follow proper surgical procedures and carefully evaluate patient suitability using standard practice and individual experience. Avoid damage to the tissue expander and use a sterile backup in case of damage. Pay careful attention to tissue tolerance and hemostasis during surgery. Expansion should proceed moderately and never beyond patient or tissue tolerance. Avoid contamination in any postoperative procedure.

ADVERSE REACTIONS

Possible adverse reactions include deflation, tissue damage, infection, extrusion, hematoma/seroma, capsular contracture, premature explantation, displacement, effects on bone, pain, sensation, distortion, inadequate tissue flap, and inflammatory reaction.

For more information, please see the full Directions for Use at www.allergan.com/products.

To report a problem with Natrelle® 133S Smooth Tissue Expanders, please call Allergan® at 1-800-624-4261.

Natrelle® 133S Smooth Tissue Expanders are restricted to sale by or on the order of a licensed physician.

Natrelle INSPIRA® Style SCX Smooth, Extra-Full Profile			
WIDTH (CM)	PROJECTION (CM)	VOLUME (CC)	Smooth Cohesive (SCX)
9.50	4.2	200	SCX-200
10.00	4.5	225	SCX-225
10.25	4.8	255	SCX-255
10.50	5.0	285	SCX-285
10.75	5.2	310	SCX-310
11.00	5.4	340	SCX-340
11.25	5.6	375	SCX-375
---	---	---	---
11.50	5.8	400	SCX-400
11.75	5.9	420	SCX-420
12.00	6.0	445	SCX-445
12.25	6.1	470	SCX-470
12.50	6.1	495	SCX-495
12.75	6.2	525	SCX-525
13.00	6.2	545	SCX-545
13.25	6.3	560	SCX-560
13.50	6.4	580	SCX-580
13.75	6.4	615	SCX-615
14.00	6.5	650	SCX-650
14.25	6.6	700	SCX-700
14.50	6.7	750	SCX-750
---	---	---	---
15.00	6.7	800	SCX-800
---	---	---	---
---	---	---	---
---	---	---	---

Natrelle® Breast Implants **IMPORTANT SAFETY INFORMATION**

WARNINGS

- Breast implants are not considered lifetime devices. The longer patients have them, the greater the chance they will develop complications, some of which will require more surgery
- Breast implants have been associated with the development of a cancer of the immune system called breast implant-associated anaplastic large cell lymphoma (BIA-ALCL). This cancer occurs more commonly in patients with textured breast implants than smooth implants, although rates are not well defined. Some patients have died from BIA-ALCL
- Patients receiving breast implants have reported a variety of systemic symptoms, such as joint pain, muscle aches, confusion, chronic fatigue, autoimmune diseases, and others. Individual patient risk for developing these symptoms has not been well established. Some patients report complete resolution of symptoms when the implants are removed without replacement

Please see additional Important Safety Information on back.

Natrelle®

Allergan
Aesthetics
an AbbVie company

© 2022 AbbVie. All rights reserved. NATRELLE® and its design are trademarks of Allergan, Inc., an AbbVie company. Natrelle.com NAT64569-v6 05/22 016849